English version 

SNAP* Leishmania

The Canine Leishmania Antibody Test Kit is an Enzyme-linked Immunosorbent assay for in vitro detection of *Leishmania infantum* (*L. infantum*) antibodies in canine whole blood, serum, or plasma.

Precautions and warnings

- All waste should be properly decontaminated prior to disposal.
- Do not mix components from kits with different serial numbers.
- Do not use a SNAP device that has been activated prior to the addition of sample.
- Refer to country specific Safety Data Sheet for regional hazard identification.

Storage

- Store at 2–8°C until expiration date.
- SNAP devices and reagents can be stored at room temperature (18–25°C) for 6 months or until the expiration date, whichever occurs first.
- After SNAP devices and reagents are removed from 2–8°C for **more than 24 hours**, the expiration date is 6 months or the printed expiration date, whichever occurs first. If the 6 month expiration date occurs prior to the printed expiration date, record the new date in the space provided on the kit.

Kit components

Item	Reagent	Quantity
1	Leishmania Antigen:HRPO Conjugate	1 x 7,0 mL
2	SNAP Device	10 or 30
Reagents contained in each device:		
Wash solution		0,4 mL
Substrate solution		0,6 mL
Other components: transfer pipettes, sample tubes and reagent rack		

Conjugate - H317/H412/P261/P272/P273/P280/P333 + P313/P362 + P364/P501: May cause an allergic skin reaction. Harmful to aquatic life with long lasting effects. Avoid breathing mist/vapours. Contaminated work clothing should not be allowed out of the workplace. Avoid release to the environment. Wear protective gloves. If skin irritation or rash occurs: Get medical advice/attention. Take off contaminated clothing and wash it before reuse. Dispose of contents/container in accordance with local/regional/national/international regulations.




WARNING  **Wash Solution** - H317/H319/H412/P261/P273/P280/P305 + P351 + P338/P333 + P313/P362 + P364: May cause an allergic skin reaction. Causes serious eye irritation. Harmful to aquatic life with long lasting effects. Avoid breathing mist/vapours. Avoid release to the environment. Wear eye protection/ face protection. Wear protective gloves. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If skin irritation or rash occurs: Get medical advice/attention. Take off contaminated clothing and wash it before reuse.

Sample information

- Samples must be at room temperature (18–25°C) before beginning test procedure.
- Serum, plasma, or anti-coagulated whole blood (e.g., EDTA, heparin), either fresh or stored at 2–8°C for up to one week, can be used.
- For longer storage, serum, or plasma can be frozen (-20°C or colder) and then recentrifuged before use.
- Hemolyzed or lypemic samples will not affect results.

Test procedure

- If stored in a refrigerator, allow all components to equilibrate at room temperature (18–25°C) for 30 minutes before use. **Do not heat.**

Positive control

Leishmania Sample Spot

Negative result

Only positive control spot develops color.

Positive result

Positive control spot and sample spot develop color.

Invalid results

- Background**
If the sample is allowed to flow past the Activate Circle, background color may result. Some background color is normal. However, if colored background obscures test result, repeat the test.
- No color development**
If positive control does not develop color, repeat the test.

Interpreting test result

To determine test results, read the reaction spots in the Result Window. Color development in Leishmania sample spot is proportional to the Leishmania-specific antibody concentration in the sample. If no color develops in the positive control spot, repeat the test.

Positive control

Leishmania Sample Spot

Negative result

Only positive control spot develops color.

Positive result

Positive control spot and sample spot develop color.

Invalid results

- Background**
If the sample is allowed to flow past the Activate Circle, background color may result. Some background color is normal. However, if colored background obscures test result, repeat the test.

- No color development**
If positive control does not develop color, repeat the test.

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

Symbol descriptions

	Use by date		Authorized representative in the European Community
	Batch code (lot)		Authorized representative in the European Community
	Serial number		Consult instructions for use
	Temperature limitation		In vitro diagnostic
	Manufacturer		Date of manufacture
	Catalog number		

IDEXX Technical Support
USA/Canada: 1 800 248 2483 • idexx.com
Australia: 1300 44 33 99 • idexx.com.au
Europe: idexx.eu

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© 2024 IDEXX Laboratories, Inc. Alle rechten voorbehouden.
© 2024 IDEXX Laboratories, Inc. Všechny práva vyhrazena.