SNAP® Combo Plus

SNAP Combo Plus FeLV Ag/FIV Ab test is an enzyme-linked immunosorbent assay (ELISA) for the simultaneous detection of feline leukemia virus (FeLV) antigen and antibodies control the virus de l'immunodéficience féline (FIV) in feline serum, plasma or whole blood.


Precautions and warnings.

- All waste should be properly decontaminated prior to disposal.
- Do not mix components from kits with different serial numbers.
- Do not use a SNAP device that has been activated prior to the addition of sample.
- Refer to country specific Safety Data Sheet for regional hazard identification.

Storage.

- SNAP devices and test reagents must be stored at 2-8°C.
- All components must be at room temperature 18-25°C before running the test - **DO NOT HEAT.**

Kit components				
Item	Reagent	5 test Quantity	15 test Quantity	30 test Quantity
1	1 bottle Anti-FeLV/FIV Ag-HRPo conjugate	3.5 mL	3.5 mL	7.0 mL
2	SNAP device: Monoclonal antibodies to p27, inactivated FIV antigen, and positive and negative controls	5	15	30
Reagents contained in each device:				
Wash solution		0.4 mL	0.4 mL	0.4 mL
Substrate solution		0.6 mL	0.6 mL	0.6 mL
Other components: transfer pipettes, sample tubes, and reagent rack				


Sample information

- Samples must be at room temperature 18-25°C before beginning test procedure.
- Serum or plasma, either fresh, previously frozen or stored at 2-8°C, may be used in this test.
- Serum or plasma may be stored for up to 7 days at 2-8°C. For longer storage, sample should be frozen (-20°C or colder).
- Previously frozen or older samples must be centrifuged before use.
- Hemolyzed or lipemic samples will not affect results.
- EDTA or heparin in plasma will not affect results.
- Whole blood may be used. Whole blood must be anticoagulated (eg. EDTA, heparin) and may be used either fresh or after refrigeration 2-8°C for up to one week.

Test procedure

- Using the pipette that is provided, transfer **3 drops of sample** (whole blood, serum, or plasma) into sample tube.
- Holding bottle vertical, add **4 drops of conjugate** to the sample tube.
- Cap the sample tube and mix thoroughly by **inverting tube 3-5 times**.
- Place the device on a flat surface. Add contents of sample tube to Sample Well, being careful not to splash contents outside of Sample Well. Sample will flow across Result Window, reaching Activate Circle in approximately 30-60 seconds. Some sample may remain in Sample Well.
Watch carefully for sample or blue color in the Activate Circle.
- When color **FIRST** appears in Activate Circle push Activator firmly until it is flush with the device body.

Note: Some samples may not flow to the activate circle within 60 seconds, and, therefore, the color may not turn color. In this case, press the Activator after sample has flowed across Result Window.

- Keep the device horizontal to ensure accurate results.
- Read test result at **10 minutes**. Note: positive control may develop sooner, but results are not complete until 10 minutes.

Interpreting test result

To determine test results, read the reaction spots in the Result Window. Color development in sample spots is proportional to the concentration of FeLV antigen or FIV antibody in the sample. If no color develops in the positive control spot, repeat the test.

Positive result	Invalid result
Only positive control spot develops color.	
Positive result	Invalid result
FeLV Antigen	FeLV antigen and FIV antibody
Positive control spot and sample spot develop color.	Positive control spot and sample spot develop color.
Positive control spot and sample spot develop color.	Positive control spot and sample spot develop color.

Reaction with negative control

The negative control spot serves as a safeguard against false positives and helps indicate that the assay has been run properly.

Positive result	Invalid result
If color in the FIV Ab or FeLV Ag sample spot is darker than negative control spot, result is positive for that spot.	If color in the negative control spot is equal to or darker than FIV Ab or FeLV Ag sample spot, result is invalid for that sample spot.

- Background**
If the sample is allowed to flow past the Activate Circle, background color may result. Some background color is normal. However, if colored background obscures test result, repeat the test.
- No color development**
If positive control does not develop color, repeat the test.

IDEXX SNAPshot Dx® Analyzer


Test results can also be read using the SNAPshot Dx Analyzer. A complete description of how to enter patient data and read test results using the SNAPshot Dx Analyzer can be found in the SNAPshot Dx Analyzer Operator's Guide.

Symbol descriptions	Descriptions des symboles
LOT Use by date	ECOREP Autorized representative in the European Community
SN Batch code (lot)	U Consult instructions for use
IVD Temperature limitation	IVD In-vitro diagnostic
M Manufacturer	D Date of manufacture
REF Catalog number	D Date of manufacture

Assistance technique IDEXX
E.U./Canada: 1 800 248 2483 • idexx.com
Australia: 1300 44 33 99 • idexx.com.au
Europe: idexx.eu

Assistenza tecnica IDEXX
E.U./Canada: 1 800 248 2483 • idexx.com
Australia: 1300 44 33 99 • idexx.com.au
Europe: idexx.eu

Assistenza tecnica IDEXX
E.U./Canada: 1 800 248 2483 • idexx.com
Australia: 1300 44 33 99 • idexx.com.au
Europe: idexx.eu

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